



## Inter-Therm™ Pro HMEF



**Airway Management** ▪ Heat and moisture exchanging filters (HMEFs)

## Inter-Therm™ Pro HMEF

The Inter-Therm Pro HMEF is a high efficiency HEPA filter combined with a unique corrugated paper HME which maximises the surface area to help provide a high level of protection combined with excellent humidification.

The Inter-Therm Pro HMEF pleated mechanical filter provides bacterial and viral efficiency of  $\geq 99.9999\%$ , reducing the risk of cross contamination, combined with optimal humidification of gases during mechanical ventilation.

### Tethered luer cap

For safe monitoring of respiratory gases

### Pleated mechanical filter

Provides high efficiency HEPA filtration

### Corrugated paper HME

Provides optimal humidification of respiratory gases

### Colour coding

Enhances identification, ensuring that the correct filter/HME is used for the correct clinical application



<b>Code</b>	<b>1791902</b>
Box Qty.	50
Moisture loss	8mg/L H <sub>2</sub> O/L
Calculate moisture return	30.6mg H <sub>2</sub> O/L
Filtration efficiency	99.9999%
Resistance at 30L/min	1.3cm H <sub>2</sub> O
Resistance at 60L/min	2.5cm H <sub>2</sub> O
Compressible volume	125ml
Weight	49g
Connectors	22M/15F-22F/15M
Minimum tidal volume	$\geq 300$ ml

Make an enquiry

IS12.28 • Issue 2 11.23



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK  
T: +44 (0)118 965 6300 info@intersurgical.com www.intersurgical.com



The manufacturer Intersurgical Ltd is certified to ISO 14001:2015, ISO 9001:2015, ISO 13485:2016 and MDSAP

**Please think before you print**  
Save energy and paper.  
If you must print this information sheet please print it double sided.