



Laryngoscope handles and blades

Combining safety and strength




Airway Management ▪ Airway Devices



Reprocessable laryngoscope handles

As part of our expanding range of airway devices and accessories, Intersurgical has introduced a range of reprocessable fibreoptic laryngoscope handles. All handles are ISO 7376-3 compliant, providing superior illumination to incandescent bulbs and a reduced risk of mucosal burns. Produced from high quality nickle plated brass, the entire range has been designed to be reprocessed by autoclave up to 400 times.¹

Our Hyper LED range offers superior light intensity in three different profiles - standard, slim and short.



Simple and secure fitting to all ISO7376-3 specification laryngoscope handles



Hyper LED handles

Standard handle
Most popular choice for routine cases

Short handle
For cases where the chest and / or breasts come into contact with handle during use, or for when cricoid pressure is being applied²

Slim handle
Product of choice for paediatrics

Code	Description	Box Qty.
7111000	Hyper LED standard reprocessable laryngoscope handle	1
7112000	Hyper LED short reprocessable laryngoscope handle	1
7113000	Hyper LED slim reprocessable laryngoscope handle	1

[Make an enquiry](#)

References:

1. In house testing.

2. Dorsch J.A. & Dorsch, S.E. 2012. Understanding Anesthesia Equipment. 5th edition. Lippincott Williams & Wilkins. Philadelphia.



Metal single use laryngoscope blades

A comprehensive range of single-use, stainless steel laryngoscope blades in Macintosh and Miller profiles providing the strength of a reusable metal blade, without the risk of cross-infection.

All Intersurgical metal blades are individually packaged and are compatible with ISO7376-3 green specification handles, as indicated by the green heel of the blade.



Macintosh style

Code	Description	Size	Box Qty.
7145000	Metal laryngoscope blade, fibre optic, Macintosh	5	5
7144000	Metal laryngoscope blade, fibre optic, Macintosh	4	40
7143000	Metal laryngoscope blade, fibre optic, Macintosh	3	40
7142000	Metal laryngoscope blade, fibre optic, Macintosh	2	10
7141000	Metal laryngoscope blade, fibre optic, Macintosh	1	10

Make an enquiry

Miller style

Code	Description	Size	Box Qty.
7242000	Metal laryngoscope blade, fibre optic, Miller	2	40
7241000	Metal laryngoscope blade, fibre optic, Miller	1	10
7240000	Metal laryngoscope blade, fibre optic, Miller	0	10


Make an enquiry





Plastic single use laryngoscope blades

Intersurgical plastic single use laryngoscope blades provide a safe, strong, light and cost-effective alternative to reusable metal versions. The range of Miller and Macintosh blades provides effective illumination and visualisation of the airway to aid with intubation and airway clearance. Our blades are compatible with the ISO7376-3 green specification handles, as indicated by the green dot near the heel of the blade.



High optical light transmission

Provides the user with an excellent light transfer to the blade tip

Single construction blade made from high impact grade ABS plastic

Rigid blade which is a cost effective, disposable alternative to reusable and disposable metal blades

Engagement to handle click mechanism via two plastic lugs

Provides secure and positive connection of the blade to the handle



Ultra smooth blade with matt finish

Smooth transition into the airway and easy passage of the tube, with no light reflection off the blade

Macintosh style



Code	Description	Size	Box Qty.
7044000	Plastic laryngoscope blade, fibre optic, Macintosh	4	40
7043000	Plastic laryngoscope blade, fibre optic, Macintosh	3	40
7042000	Plastic laryngoscope blade, fibre optic, Macintosh	2	50

Miller style



Code	Description	Size	Box Qty.
7041000	Plastic laryngoscope blade, fibre optic, Miller	1	70
7040000	Plastic laryngoscope blade, fibre optic, Miller	0	70

Make an enquiry



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK

T: +44 (0)118 965 6300 F: +44 (0)118 965 6356 info@intersurgical.com www.intersurgical.com



The manufacturer Intersurgical Ltd is certified to ISO 9001:2015, ISO 13485:2016 and ISO 14001:2015

Please think before you print
Save energy and paper.
If you must print this information sheet please print it double sided.

IS6.11 • Issue 3 07.19