


# Heated Wire Breathing Systems

for use in Adult Intensive Care

22mm Heated Wire Breathing Systems



22mm Heated Wire Breathing Systems




Intersurgical are delighted to launch a new range of heated wire breathing systems to add to the existing range of Flextube and Smoothbore products. A new connector has been developed to allow connection directly to the Fisher and Paykel MR850 humidifier base, these systems can also be used on the Fisher and Paykel MR730 humidifier bases if required. All the heated wire systems are available with single or dual heated wire limbs, three in one option or CPAP alternative.


Intersurgical's extensive research in many different ITU conditions has proven that our new compound wire combined with the insulating sleeve on the inspiratory limb provides optimum performance and reduces the 'rain out' formed within the system to an absolute minimum.

22mm Standard Breathing Systems

Code: 2025	Single Heated Wire System 1.6m length + 0.8m limb	Length: 1.6m	Box Qty: 10
------------	---	--------------	-------------




Code: 2019	Breathing System - 1.6m + water trap	Length: 1.6m	Box Qty: 10
------------	--------------------------------------	--------------	-------------



Smoothbore Heated Wire Breathing Systems


Code: 5025   1.6m length single heated wire system with water trap

Box Qty: 10




Code: 5019   1.6m length single heated wire system with water trap

Box Qty: 10



Code: 5026   1.6m length dual heated wire system


Box Qty: 10



22mm Breathing Systems for use with Fisher & Paykel Humidifiers


Code: 2026 Dual Heated Wire System 1.6m length + 0.8m limb

Box Qty: 10




Code: 2039 Dual Heated Wire System 1.6m length + 0.8m limb

Box Qty: 10




Adult 22mm Flextube™ Breathing System for Ventilation → CPAP → Aerosol Therapy


Code: 2027	MR850 Option 1.6m length + 0.8 limb (product shown below)	Box Qty: 10
Code: 2049	1.6m length + 0.8 limb	Box Qty: 15



Or CPAP...




Or Aerosol Therapy





CPAP Breathing Systems with Heated Wire

Code: 2050	MR850 Option 2.4m length + 1.2m limb (product shown below)	Length: 1.6m	Box Qty: 10
Code: 2037	22mm Dräger Breathing System	Length: 1.6m	Box Qty: 20



Code: 2038	1.75m length + 1.6m limb	Length: 1.6m	Box Qty: 20
------------	--------------------------	--------------	-------------





## Intersurgical Offer:

A Wide Product Range  
High Quality Products  
All Products CE Marked  
Fast, Efficient Service  
Reliability  
Excellent Value

### Deutschland

T 02241 311063 F 02241 313143  
[info@intersurgical.de](mailto:info@intersurgical.de)

### France

T 01 48 76 72 30 F 01 48 76 91 38  
[info@intersurgical.fr](mailto:info@intersurgical.fr)

### España

T 91 665 73 15 F 91 685 50 99  
[info@intersurgical-es.com](mailto:info@intersurgical-es.com)

### Portugal

T 219108550 F 219108559  
[info@intersurgical2.pt](mailto:info@intersurgical2.pt)

### Lithuania

T 37017 66611 F 37017 66622  
[info@intersurgical.lt](mailto:info@intersurgical.lt)

### Nederland

T 0413 243860 F 0413 243869  
[info@intersurgical.demon.nl](mailto:info@intersurgical.demon.nl)

### United Kingdom

Crane House, Molly Millars Lane  
Wokingham RG41 2RZ

T + 44 (0)118 9656 300  
F + 44 (0)118 9656 356  
[info@intersurgical.co.uk](mailto:info@intersurgical.co.uk)  
[www.intersurgical.com](http://www.intersurgical.com)