

# Pocket Resuscitation Face Mask

The quick and convenient,  
ready to use resuscitation mask



**INTERSURGICAL**<sup>®</sup>  
COMPLETE RESPIRATORY SYSTEMS



# Pocket Resuscitation Face Mask



## The quick and convenient, ready to use resuscitation mask

When an unconscious patient is not breathing, immediate action is required. Every second counts and the sooner resuscitation is commenced the better the chance of a successful outcome.

Mouth to mouth has been the traditional means of giving immediate respiratory support. However, due to the ever increasing incidence of infectious disease and the inherent unhygienic nature of the practice there may be a natural reluctance to carry it out.

Intersurgical offer a convenient ready to use mouth to mask device that allows you to give immediate respiratory support to your patient without compromising yourself.

The easily identifiable protective case allows the Pocket Resuscitation Face Mask to be stored in an accessible area in the department,

included in an emergency pack or to be carried for personal use.

The crystal clear mask allows clear visualisation of the patient's colour giving an immediate indication of the effectiveness of the respiratory support. Any vomit, mucous or foreign body blocking the patient's airway can be easily seen. The soft cuff moulds itself to the contours of the patient's face providing a comfortable air tight seal.


The one way protection valve allows gas to pass from the resuscitator to patient without any physical contact between the two. This will minimise any risk of cross contamination and vastly improves the hygiene of mouth to mouth resuscitation.

The mouthpiece is a standard 15mm male connection so it will fit our standard resuscitation systems and breathing systems.

Code: 1157	Pocket Resuscitation Face Mask Replacement Valve	Box Qty: 5
------------	--	------------



Code: 1158	Pocket Resuscitation Face Mask with supplementary oxygen connector, headstrap and filter	Box Qty: 16
------------	--	-------------



For details of our full range of Anaesthetic Breathing Systems and Accessories please refer to our full product catalogue or the Intersurgical website

## Intersurgical Offer:

- A Wide Product Range
- High Quality Products
- All Products CE Marked
- Fast, Efficient Service
- Reliability
- Excellent Value

### Deutschland

T 02241 311063 F 02241 313143  
info@intersurgical.de

### France

T 01 48 76 72 30 F 01 48 76 91 38  
intersurgical@wanadoo.fr

### España

T 91 653 96 11 F 91 653 99 23  
intersurg-es@nexo.es

### Portugal

T 01 498 98 80 F 01 498 98 89  
mop29707@mail.telepac.pt

### Lithuania

T 37017 54541 F 37088 52773  
enquiries@intersurgical.lt

### Nederland

T 0413 243860 F 0413 243869  
info@intersurgical.demon.nl

### United Kingdom

Crane House, Molly Millars Lane  
Wokingham RG41 2RZ

T + 44 (0)118 9656 300  
F + 44 (0)118 9656 356  
enquiries@intersurgical.co.uk  
www.intersurgical.com