

One-piece® Guedel airways

Safer by design




Safer by Design

Utilisation of the latest manufacturing technology has enabled us to produce the first major breakthrough in Guedel Airway design for thirty years.

We have combined two materials with very different physical characteristics into one moulding. The result is the elimination of the separate bite block found in all conventional Guedel Airway designs.

Loose or detached bite blocks constitute a real threat to patient safety and have been the subject of several hazard notices in the past. Our One-piece design overcomes all of the associated risks.

As well as improved safety, the One-piece also incorporates a soft tip reducing the risk of trauma to the patient.



Ordering Information

1115120 (S*)	Guedel Airway size 5 (12.0), purple colour-coding	Box qty. 70
1114100 (S*)	Guedel Airway size 4 (10.0), red colour-coding	Box qty. 80
1113090 (S*)	Guedel Airway size 3 (9.0), orange colour-coding	Box qty. 80
1112080 (S*)	Guedel Airway size 2 (8.0), green colour-coding	Box qty. 100
1111570 (S*)	Guedel Airway size 1.5 (7.0), yellow colour-coding	Box qty. 100
1111065 (S*)	Guedel Airway size 1 (6.5), white colour-coding	Box qty. 20
1110055 (S*)	Guedel Airway size 0 (5.5), grey colour-coding	Box qty. 25
1100050 (S*)	Guedel Airway size 00 (5.0), blue colour-coding	Box qty. 25
1000035 (S*)	Guedel Airway size 000 (3.5), pink colour-coding	Box qty. 25

* Add an S to the code number for the sterile version of this product.

IS10.15 Issue 4 MD 09.12